

Protocol 3H - Waveform Capnography Adult & Pediatric

Changes were made to under Technique:

Physio-Control LifePak[®] 12/15) to Monitor EtCO₂:

1. Make sure the monitor is ON.
2. Select the appropriate EtCO₂ accessory for the patient.
3. Open the CO₂ port door and insert the FilterLine[®] connector; turn connector clockwise until tight.
4. Verify that the CO₂ area is displayed. The EtCO₂ monitor performs the autozero routine as part of the initialization self-test.
5. Display CO₂ waveform in Channel 3 on the LifePak screen.
6. Connect the CO₂ FilterLine[®] set to the patient.
7. Confirm that the EtCO₂ value and waveform are displayed.

NOTE: Do not connect the FilterLine[®] set to the patient/ventilation system until the EtCO₂ monitor has completed its self-test and warm-up.

Protocol 5A: Chest Pain - Uncertain Etiology Adult & Pediatric

Changes to the current treatment protocol:

Added to the Treatment Priorities under 5 in 10 minutes of patient contact number:

1. Adult- ≥ 35 years of age - ASA ; in the EMR/EMT treatment box it was also added:
Adult: ≥ 35 years of age - ASA 324/325 mg.

Protocol 16E – Aspirin

Changes to the current protocol:

Added to the Indications: Chest Pain - Uncertain Etiology (5A) ≥ 35 years of age. This was added to support the treatment in Protocol 5A.