


4N – RESUSCITATION TEAM ROLES EleGARD® HEAD UP CPR-ADULT

TREATMENT PRIORITIES PRIOR TO PRIMING PHASE < 2 MINUTES:

1. Circulatory Support
 - Apply ResQCPR®, and place EleGARD®
 - Continuous chest compression: 110/min standard CPR
 - Place i-gel® with ResQPOD®, inline EtCO2/waveform capnography, & OGT
 - Once i-gel®, and ResQPod® in place: 80 chest compression – decompression cycles/min ResQCPR® with ResQPUMP®
 - Prior to priming phase, rhythm analysis, and timely defibrillation (if indicated)
 - Start priming phase by pressing button # 2
 - Priming phase is 2 min uninterrupted CPR with, i-gel®, ResQPod®, ResQPUMP®, and EleGARD®
2. Oxygenation/Ventilation support
 - Avoid hyperventilation in rate & volume
 - Use waveform capnography **mandatory use if patient intubated



Post Deployment of EleGARD® Follow

- Protocol 4B - Resuscitation Team Roles
- 4N – EleGARD® HEAD UP CPR ADULT